





The real-time setup

Michael Lührs





Overview

1. General real-time fMRI setup
2. Potential pitfalls
3. Important considerations



1. General real-time fMRI setup






Individual real-time fMRI components

- MRI
 - Siemens
 - Phillips
 - GE
 - United Imaging
- Real-time fMRI processing computer
 - Windows, Linux, Mac
- Stimulus/Neurofeedback presentation computer
 - Windows, Linux, Mac
- Peripherals
 - Physiological recordings, Motion, Eye Tracking, Response box, ...



Real-Time fMRI setup





Real-Time fMRI setup

MRI console computer




IP: 192.168.1.1

Real-time fMRI computer



IP: 192.168.1.2

Network connection



Both computers have to be in the same network to access the shared folder provided by the Turbo-BrainVoyager computer. The folder can be reached on the defined location.


Shared folder "rtfmri"



Location: \\192.168.1.2\rtfmri



Neurofeedback setup





Providing feedback





TCP connection




- On Stimulation PC:
 - TBV Network Interface

- On Analysis PC:
 - TBV Network Plugin



2. Potential pitfalls





Real-Time fMRI setup

MRI console computer




IP: 192.168.1.1

Real-time fMRI computer



IP: 192.168.1.2

Network connection



Both computers have to be in the same network to access the shared folder provided by the Turbo-BrainVoyager computer. The folder can be reached on the defined location.


Shared folder "rtfmri"



Location: \\192.168.1.2\rtfmri




Network structure









Data transfer times



Lührs, Michael, et al. "Fast Retrieval of fMRI Data for Real-Time Applications: Improving the Transfer Time through Direct Connection." *Aperture Neuro*, vol. 3, Aug. 2023, <https://doi.org/10.52294/001c.77768>.




Data transfer times



Lührs, Michael, et al. "Fast Retrieval of fMRI Data for Real-Time Applications: Improving the Transfer Time through Direct Connection." *Aperture Neuro*, vol. 3, Aug. 2023, <https://doi.org/10.52294/001c.77768>.



Data transfer times



Renz, Malika P et al. "Practical challenges of continuous real-time functional magnetic resonance imaging neurofeedback with multiband accelerated echo-planar imaging and short repetition times." Human brain mapping vol. 44,3 (2023): 1278-1282. doi:10.1002/hbm.26154







Continuous vs. intermittent feedback

- Continuous:
 - Give feedback as soon as it is available
 - During regulation attempts
- Intermittent:
 - Give feedback only after regulation attempts



Real-time pre-processing and analysis

- Delays in data processing/analysis
 - Variability due to:
 - Motion parameters
 - GLM statistics
 - MVPA
 - Connectivity
 - Real-time processing computer is not made for real real-time
 - Delays can occur more in later runs compared to initial




3.

Important considerations



Verify real-time export parameters of your study

- Acquire phantom scans at least 20% longer than your planned
 - Use the same sequence parameters as for the final study
 - Record export times and repeat this test on different days and times
- Include longer baseline phase in the beginning of each run
 - Accounting for stable data reconstruction (Multi-band)
 - More stable statistics for real-time GLM
- Log transfer time during experiment
 - Important to be able to trace potential delays also month after the study



Questions?

